

**Tube for
PF6 / PD12 / FP50 / FPC50**




Sales manual

Rev.: 05

25.02.2011

AEF



Part number	Description	A (ID) (mm)	B (wall thickness) (mm)	C (length) (m)
84-103-005	ACCUSIL™ Platinum cured silicone tubing	0.5	1.6	10
84-103-008	ACCUSIL™ Platinum cured silicone tubing	0.8	1.6	10
84-103-012	ACCUSIL™ Platinum cured silicone tubing	1.2	1.6	10
84-103-016	ACCUSIL™ Platinum cured silicone tubing	1.6	1.6	10
84-103-032	ACCUSIL™ Platinum cured silicone tubing	3.2	1.8	10
84-103-048	ACCUSIL™ Platinum cured silicone tubing	4.8	2.0	10
84-103-060	ACCUSIL™ Platinum cured silicone tubing	6.0	2.1	10
84-103-080	ACCUSIL™ Platinum cured silicone tubing	8.0	2.2	10
84-105-008	Pharmed tubing	0.8	1.6	1
84-105-016	Pharmed tubing	1.6	1.6	1
84-105-032	Pharmed tubing (32-152-100)	3.2	1.6	1
84-105-048	Pharmed tubing (32-152-101)	4.8	1.6	1
84-105-060	Pharmed tubing (32-152-102)	6.0	1.6	1
84-105-080	Pharmed tubing (32-152-103)	8.0	1.6	1
84-301-005	Silicone tubing	0.5	1.6	5
84-301-008	Silicone tubing	0.8	1.6	5
84-301-012	Silicone tubing	1.2	1.6	5
84-301-016	Silicone tubing	1.6	1.6	5
84-301-032	Silicone tubing	3.2	1.6	5
84-301-048	Silicone tubing	4.8	2.0	5
84-301-060	Silicone tubing	6.0	2.0	5
84-301-080	Silicone tubing	8.0	2.0	5
84-301-200	Set of silicone tubing from 84-301-005 to 84-301-080			
84-311-008	Silicone tubing	0.8	1.6	25
84-311-016	Silicone tubing	1.6	1.6	25
84-311-032	Silicone tubing	3.2	1.6	25
84-311-048	Silicone tubing	4.8	2.0	25
84-311-060	Silicone tubing	6.0	2.0	25
84-311-080	Silicone tubing	8.0	2.0	25

- Material certificate is supplied with all tubing, except for the pharmed tubing.
- Tube locks for 3.2 – 4.8 – 6.0 and 8.0 pharmed tubing are not supplied together with the peristaltic dispenser. Part numbers for tube locks are mentioned in the brackets.