


RT60/80_CT60/80

ROTARY TABLES



- Rotary bottle feeding table
- Rotary bottle collecting table
- Two disc diameters: 600 mm or 800 mm
- Optional mounting of collecting plate on front or rear of rotary table (ie. right or left-handed bottle in- and out-feed)

Watson-Marlow Flexicon

Frejasvej 2-6
DK-4100 Ringsted
Denmark
Tel. +45 57 67 11 55
Fax. +45 57 67 05 41
flexicon@flexicon.dk
www.wmflexicon.dk




RT60/80

TECHNICAL SPECIFICATIONS

CT60/80

Rotary feeding table for vial accumulation - type RT60 or RT80



Here the collecting plate is shown as front mounted (right-handed infeed). Optional mounting of collecting plate to the rear (left-handed infeed) as well.

Disc sizes:

600 mm or 800 mm diameter

Transport directions:

Right or left-handed bottle infeed or bottle outfeed

Disc speed:

Adjustable

Vial/bottle size:

Up to 100 mm diameter and 250 mm height

Adjustment to bottle sizes:

Easily adjusted to various bottle sizes with no use of tools – for collection tables no adjustments are required

Cleaning:

Easy to clean - all parts in contact with vials/bottles are in stainless steel – AISI304

Electrical supply:

240 V single phase electrical supply - or connection to Watson-Marlow Flexicon monobloc system

Rotary collecting table for vial accumulation - type CT60 or CT80

