


DN20

OEM DIVING NOZZLE SYSTEM



- Diving nozzle support
- 1, 2 or 4 filling nozzle
- Compact and user-friendly
- Fully automatic when combined with Watson-Marlow Flexicon multi filling system or gating system

Watson-Marlow Flexicon

Frejasvej 2-6
DK-4100 Ringsted
Denmark
Tel. +45 57 67 11 55
Fax. +45 57 67 05 41
flexicon@flexicon.dk
www.wmflexicon.dk



DN20 TECHNICAL SPECIFICATIONS

The Watson-Marlow Flexicon DN20 diving nozzle system is a very compact and user-friendly filling nozzle support, which dives the nozzle into the bottle during the fill. It can be mounted onto any conveyor and can support up to 4 nozzles at a time.

It is ideal in combination with:

- Watson-Marlow Flexicon Multi Filling System
- Watson-Marlow Flexicon gating system

These combinations will create a fully automatic filling system for bottles of various sizes.

The diving system combined with the Watson-Marlow Flexicon gating system and filling units create a fully automatic filling system.

This system is able to fill both liquids like water or more viscous products like shampoo or oil. This all depends on which of the pumps you choose from the Watson-Marlow Flexicon Multi Filling System.

Number of filling nozzles:

1, 2 or 4

Weight:

16 kg



The Watson-Marlow Flexicon Filling System can be a combination of any of the above machines.